



Best Practices for Reducing the Risk of SSI's Using an Infection Prevention Bundle

Melissa Kelly, RN, BSN, CNOR
Medical Liaison

October 9, 2025

© Solventum 2024. All rights reserved.



Important Information

Prior to the use of any Solventum and/or 3M Therapy System, it is important for the provider to consult the treating physician and read and understand all Instructions for Use, including Safety Information, Dressing Application Instructions, and Therapy Device Instructions.

Specific indications, contraindications, warnings, precautions, and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

To the extent this presentation contains case studies and clinical reports, the results and outcomes should not be interpreted as a guarantee or warranty of similar results. Individual results may vary depending on the patient's circumstances and condition.

This information is intended for healthcare professionals only. Solventum recommends that clinicians participate in device in-service and training prior to use.

Solventum and/or 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.

Follow local institutional protocols for infection control and waste disposal procedures. Local protocols should be based on the applicable federal, state and/or local government environmental regulations.


Objectives

Explain the SSI Risk Equation.

Discuss Solventum solutions as part of a bundle to help reduce the dose of bacteria on a surgical patient.

Discuss a Solventum solution in the bundle that can help improve the resistance of the patient.

Surgical Site Infections (SSIs)



Inpatient surgical procedures ~10 million in U.S.¹


SSIs occur in 2 to 5% of patients undergoing inpatient surgery²⁻⁴

Majority of SSIs are considered preventable with evidence-based guidelines²⁻⁴

1. Surgeries in Hospital-Based Ambulatory Surgery and Hospital Inpatient Settings, 2014 #223. Hcup-us.ahrq.gov. <https://www.hcup-us.ahrq.gov/reports/statbriefs/sb223-Ambulatory-Inpatient-Surgeries-2014.jsp>. Revised July 2020. Accessed January 20, 2023.
2. Surgical Site Infections. Psnet.ahrq.gov. <https://psnet.ahrq.gov/primer/surgical-site-infections>. Published 2019. Accessed January 20, 2023.
3. Anderson DJ, Podgorny K, Berrios-Torres SI, et al. Strategies to prevent surgical site infections in acute care hospitals: 2014 update. *Infect Control Hosp Epidemiol*. 2014;35(6):605-627.
4. Fields AC, Pradarelli JC, Itani KMF. Preventing Surgical Site Infections: Looking Beyond the Current Guidelines. *JAMA*. 2020;323(11):1087–1088

The SSI risk equation

**Risk
of
SSI** =
$$\frac{\text{Dose of bacteria} \times \text{Virulence of bacteria}}{\text{Resistance of the host (patient)}}$$



Control what you
can control



Mangram AJ, Horan TC, Pearson ML, Silver LC, and Jarvis WR. Guideline for prevention of surgical site infection. *Infect Control Hosp Epidemiol.* 1999;4:247-278. Retrieved from <https://www.ncbi.nlm.nih.gov/pubmed/24799638>.

Reducing the load of bacteria

The SSI risk equation

Risk of SSI

$$= \frac{\begin{array}{c} \boxed{5} \\ \text{Dose of bacteria} \end{array} \times \begin{array}{c} \boxed{5} \\ \text{Virulence of bacteria} \end{array}}{\begin{array}{c} \boxed{5} \\ \text{Resistance of the host (patient)} \end{array}} = 5$$

Control what you
can control



Mangram AJ, Horan TC, Pearson ML, Silver LC, and Jarvis WR. Guideline for prevention of surgical site infection. *Infect Control Hosp Epidemiol.* 1999;4:247-278. Retrieved from <https://www.ncbi.nlm.nih.gov/pubmed/24799638>.

The SSI risk equation

Risk of SSI

$$= \frac{\begin{array}{c} \boxed{1} \\ \text{Dose of bacteria} \end{array} \times \begin{array}{c} \boxed{5} \\ \text{Virulence of bacteria} \end{array}}{\begin{array}{c} \boxed{5} \\ \text{Resistance of the host (patient)} \end{array}} = 1$$

Control what you
can control



Mangram AJ, Horan TC, Pearson ML, Silver LC, and Jarvis WR. Guideline for prevention of surgical site infection. *Infect Control Hosp Epidemiol.* 1999;4:247-278. Retrieved from <https://www.ncbi.nlm.nih.gov/pubmed/24799638>.

Evidence-based prevention strategies to reduce the risk of SSIs

Hospital Environment

- ✓ OR ventilation
- ✓ OR surfaces
- ✓ OR traffic flow
- ✓ Microbiological sampling
- ✓ Reprocessing of surgical instruments
- ✓ Sterile field management

Patient Preparation

- ✓ Antiseptic showering (bathing)
- ✓ Nasal decolonization
- ✓ Hair removal
- ✓ Antimicrobial prophylaxis
- ✓ Preoperative patient warming
- ✓ Glucose monitoring


Surgical Intervention

- ✓ Sterilization assurance
- ✓ Vascular access
- ✓ Hand hygiene
- ✓ Temperature monitoring
- ✓ Intraoperative patient warming
- ✓ Surgical skin antisepsis
- ✓ Antimicrobial incise draping

Post Operative


- ✓ Post operative incision management
- ✓ Negative pressure
- ✓ Post operative patient warming

Nasal Decolonization



Nasal colonization can increase the risk of infection

Staphylococcus aureus causes more health care associated infections (HAIs) than any other pathogen.¹



1. Weiner LM, Webb AK, Limbago B, et al. Antimicrobial-Resistant Pathogens Associated with Healthcare-Associated Infections: Summary of Data Reported to the National Healthcare Safety Network at the Centers for Disease Control and Prevention, 2011-2014. *Infect Control Hosp Epidemiol.* 2016;37(11):1288-1301.

2. Kuehnert MJ, Kruszon-Moran D, Hill HA, et al. Prevalence of *Staphylococcus aureus* nasal colonization in the United States, 2001-2002. *J Infect Dis.* 2006;193(2):172-179.

3. Perl TM, Cullen JJ, Wenzel RP, et al. Intranasal mupirocin to prevent postoperative *Staphylococcus aureus* infections. *N Engl J Med.* 2002;346(24):1871-1877.

4. Kalmeijer MD, van Nieuwland-Bollen E, Bogaers-Hofman D, de Baere GA. Nasal carriage of *Staphylococcus aureus* is a major risk factor for surgical-site infections in orthopedic surgery. *Infect Control Hosp Epidemiol.* 2000;21(5):319-323.

5. Kluytmans JA, Mouton JW, Ijzerman EP, et al. Nasal carriage of *Staphylococcus aureus* as a major risk factor for wound infections after cardiac surgery. *J Infect Dis.* 1995;171(1):216-219.

6. Wertheim HF, Vos MC, Ott A, et al. Risk and outcome of nosocomial *Staphylococcus aureus* bacteraemia in nasal carriers versus non-carriers. *Lancet.* 2004;364(9435):703-705. doi:10.1016/S0140-6736(04)16897-9

7. Anderson DJ, Kaye KS, Chen LF, et al. Clinical and financial outcomes due to methicillin resistant *Staphylococcus aureus* surgical site infection: a multi-center matched outcomes study. *PLoS ONE.* 2009;4(12):e8305.

Decolonization product selection

Challenges of the nasal environment:

- Dark
- Moist
- High levels of microorganisms^{1,2}
- Anatomical challenges – tip of the nose, nasal orifices^{1,2}
- Cilia and mucociliary clearance³
- pH range 5.5-6.5⁴

**The goal –
select an effective
formulation that can
address challenges**

1. Sakr A, Brégeon F, Mège JL, Rolain JM, Blin O. *Staphylococcus aureus* Nasal Colonization: An Update on Mechanisms, Epidemiology, Risk Factors, and Subsequent Infections. *Front Microbiol.* 2018;9:2419. Published 2018 Oct 8. doi:10.3389/fmicb.2018.02419
2. Kumpitsch C, Koskinen K, Schöpf V, Moissl-Eichinger C. The microbiome of the upper respiratory tract in health and disease. *BMC Biol.* 2019;17(1):87. Published 2019 Nov 7. doi:10.1186/s12915-019-0703-z
3. Bustamante-Marin XM, Ostrowski LE. Cilia and Mucociliary Clearance. *Cold Spring Harb Perspect Biol.* 2017;9(4):a028241. Published 2017 Apr 3. doi:10.1101/cshperspect.a028241
4. England RJ, Homer JJ, Knight LC, Ell SR. Nasal pH measurement: a reliable and repeatable parameter. *Clin Otolaryngol Allied Sci.* 1999;24(1):67-68. doi:10.1046/j.1365-2273.1999.00223.x

Product options for nasal decolonization

- Mupirocin
- 5% Povidone Iodine
- 10% Povidone Iodine
- Alcohol

Antiseptic versus antibiotic comparison

Antiseptic	Antibiotic
Biocides - destroy or inhibit the growth of microorganisms ¹	Can either induce cell death (bactericidal drugs) or merely inhibit cell growth (bacteriostatic drugs) ⁴
Help reduce risk of infection ¹	Treat or prevent infection ³
Broad spectrum ¹	Often narrower spectrum ^{1,3}
Kill bacteria by attacking multiple cell processes ^{1,2}	Kills or inhibits growth of bacteria by inactivation of a single target ²
Nontoxic at relatively high concentrations ²	Can cause toxicity at higher concentrations, especially systemically
Resistance usually does not form ¹	Resistance can form ³

1. McDonnell G. Block's Disinfection, Sterilization, and Preservation. In: Lippincott Williams & Wilkins; 2020.

2. McDonnell G, Russell AD. Antiseptics and disinfectants: activity, action, and resistance [published correction appears in Clin Microbiol Rev 2001 Jan;14(1):227]. *Clin Microbiol Rev.* 1999;12(1):147-179.

3. Kohanski MA, Dwyer DJ, Collins JJ. How antibiotics kill bacteria: from targets to networks. *Nat Rev Microbiol.* 2010;8(6):423-435. doi:10.1038/nrmicro2333

4. Overview of Antibiotics - Infections - Merck Manuals Consumer Version. <https://www.merckmanuals.com/en-ca/home/infections/antibiotics/overview-of-antibiotics>.

3M™ Skin and Nasal Antiseptic

(Povidone-Iodine Solution 5% w/w [0.5% available iodine] USP)

Patient Preoperative Skin Preparation



Foam applicator delivery to the nares

Allows application with gentle scrubbing

- Increases antimicrobial efficacy
- Dual-action (antiseptic + scrubbing)

Foam tipped applicators

- Absorbs and delivers effective amount of solution
- Optimized tip size for application to the nares
- Gentle feel



Unique chemistry of 3M™ Skin and Nasal Antiseptic

Active: Povidone-Iodine Solution 5% w/w¹

Formulated with:

Unique film forming polymer: High viscosity gel- stays in contact with the anterior nares longer for sustained activity at the site of colonization and prevents formulation from dripping during application

Unique Buffer System:


- Maintains optimal pH
- Optimal antimicrobial activity: Iodine is rapidly reduced to inactive iodide at pH > 4.5 and hence product is buffered to help protect the reduction from happening to maintain the active iodine¹
- Product safety (minimal irritation)
- Formulation stability (shelf-life)

Combination of Surfactants: Allows for maintaining physical and chemical stability of composition and for improved spreading to completely cover the anterior nares



1. Solventum Data on File

3M™ Skin and Nasal Antiseptic



Reduces bacterial counts in the nares in 1 hour, including *Staphylococcus aureus* by 99.5% and maintains this reduction for at least 12 hours.*¹



Reduces bacteria on the groin and abdomen in 10 minutes*²



3M Skin and Nasal Antiseptic can be used on children as young as 2 months of age^{‡3}



Supports antibiotic stewardship

Iodine allergies

- There is no correlation between shellfish allergies and iodine.^{1,2}
- Iodine is a trace element essential to life and present throughout the body. A true allergy to the trace element iodine does not exist.
- An anaphylactic reaction to topical iodine antiseptic solutions is exceedingly rare and not proven to be related to iodine.
- Contact dermatitis related to topically applied iodine antiseptics is not an indication of an iodine allergy; rather, it is indicative of a reaction to chemicals in the product.^{3,4}
 - Today, topically applied iodine antiseptics contain a carrier complex (e.g. povidone, povacrylex). In very rare cases, these antiseptics may have an allergic potential.
- Stated iodine allergy warrants further conversations and education

1. Iodine “allergy” and fish allergy. February 24, 2020. *American Academy of Allergy Asthma and Immunology*. <https://www.aaaai.org/allergist-resources/ask-the-expert/answers/old-ask-the-experts/iodine-allergy-fish>.
2. Shellfish allergy is not a shell game. January 10, 2024. *American Academy of Allergy Asthma and Immunology*. <https://www.aaaai.org/tools-for-the-public/conditions-library/allergies/shellfish-allergy-can-be-dangerous>.
3. Quatresooz P, Xhaufaire-Uhoda E, Piérard-Franchimont C, Piérard GE. Regional variability in stratum corneum reactivity to antiseptic formulations. *Contact Dermatitis*. 2007;56(5):271-273. doi:10.1111/j.1600-0536.2007.01097.x
4. Sicherer SH. Risk of severe allergic reactions from the use of potassium iodide for radiation emergencies. *J Allergy Clin Immunol*. 2004;114(6):1395-1397. doi:10.1016/j.jaci.2004.09.026

CDC Strategies to Prevent Hospital-onset *Staphylococcus aureus* Bloodstream Infections in Acute Care Facilities and SSI Prevention

Patient Type	Intensive Care Unit	Non-Intensive Care Unit
Central venous catheter (CVC) or Midline Catheter Present	Topical chlorhexidine gluconate (at least 2%) + Intranasal antistaphylococcal antibiotic/antiseptic (e.g. mupirocin or iodophor) (core strategy)	Topical chlorhexidine gluconate (at least 2%) + Intranasal antistaphylococcal antibiotic/antiseptic (e.g. mupirocin or iodophor) (supplemental strategy)
No CVC or Midline Catheter present	Topical chlorhexidine gluconate (at least 2%) + Intranasal antistaphylococcal antibiotic/antiseptic (i.e. mupirocin or iodophor) (core strategy)	None (note that decolonization or pathogen reduction strategies may apply to pre-operative surgical patients outside the intensive care unit)
Surgical site infection (SSI) prevention practices	For patients undergoing high risk surgeries (e.g. cardiothoracic, orthopedic, and neurosurgery), use an intranasal antistaphylococcal antibiotic/ antiseptic (e.g. mupirocin or iodophor) and chlorhexidine wash or wipes prior to surgery (core strategy)* * Facilities can choose to apply the selected pre-operative decolonization or pathogen reduction regimen universally to all patients or can screen patients undergoing a high-risk surgery with a test that detects both MSSA and MRSA and provide the decolonization regimen only to those from whom <i>S. aureus</i> is identified.	

SHEA/IDSA/APIC

Strategies to prevent surgical site infections in acute-care hospitals: 2022 update

Risk Factor	Recommendation	Quality of Evidence
<i>S. Aureus</i> nasal colonization	Decolonize patients with nasal mupirocin or povidone-iodine prior to surgery	Moderate

Decolonization of surgical patients with an antistaphylococcal agent for cardiothoracic and orthopedic procedures was reclassified from an “additional approach” to an “essential practice”

Calderwood MS, Anderson DJ, Bratzler DW, et al. Strategies to prevent surgical site infections in acute-care hospitals: 2022 Update. Infect Control Hosp Epidemiol. 2023;44(5):695-720. doi:10.1017/ice.2023.67.

ERAS Cardiac Society

Turnkey Order Set for Surgical-Site Infection Prevention: Proceedings from the American Association for Thoracic Surgery ERAS Conclave 2023

Preoperative recommendations include screening for and treating *Staphylococcus aureus* nasal carriage.

Preoperative

- Obtain nasal swab,
 - If time permitting to get results:
 - If positive for *S. aureus*, Mupirocin 2% nasal ointment 2x/day for 5 days
 - If time not permitting to get results:
 - Empirically treat: Mupirocin 2% nasal ointment 2x/day for 5 days, discontinue if results come back negative

Or

- *Administer povidone-iodine solution 5%(0.5% available iodine) USP nasal antiseptic single application in each naris on morning of surgery*


Association of periOperative Nurses (AORN)

- Use a risk-based approach to determine *S. aureus* decolonization strategies
- Evaluate resources and plan for expansion of services for the strategy
- Select horizontal, vertical, or blended strategy for decolonization
- Determine screening method
- Determine decolonization protocol
- Implement the decolonization protocol
- Monitor the regimen

Nasal Decolonization Options Listed

- Mupirocin (22 studies)
- Povidone Iodine (7 studies)
- Octenidine (1 study)
- Alcohol based (2 studies)

Skin Prep




Goals of Antiseptic Skin Prep Solutions¹




Significantly decrease the number of bacteria colonizing on the skin, thereby reducing the risk of contamination of the incisional wound.



Antiseptics mechanically remove and chemically kill contaminating and colonizing flora.



Must provide broad-spectrum coverage.



1. Dumville JC, McFarlane E, Edwards P, Lipp A, Holmes A, Liu Z. Preoperative skin antiseptics for preventing surgical wound infections after clean surgery. Cochrane Database Syst Rev. 2015 Apr 21;2015(4):CD003949. doi: 10.1002/14651858.CD003949.pub4.

Global Clinical Guidelines

Surgical Skin Preparation



Clinical guidelines recommend to use a skin prep that contains **both** alcohol and a second active ingredient unless contraindicated

- Alcohol should not be used on mucous membranes or open wounds
- Chlorhexidine Gluconate (CHG) should not be used on mucous membranes, meninges, eyes, ears
- Allergy to CHG or Iodine

WHO	World Health Organization, 2018
AORN	Association of periOperative Registered Nurses, 2023
SHEA/IDSA	The Society for Healthcare Epidemiology of America/The Infectious Diseases Society of America, 2022
NICE	The National Institute for Health and Care Excellence, 2019
KRINKO	The Commission for Hospital Hygiene and Infection Prevention at Robert Koch Institute (2018)
CDC	Centers for Disease Control and Prevention, 2017
ACOG	American College of Obstetricians and Gynecologists, 2018
ACS/SIS	American College of Surgeons and Surgical Infection Society, 2016

Global Clinical Guidelines

Single-Use Container vs Multi-Use Bottle



Clinical guidelines recommend use of a single-use container to reduce potential contamination of skin antiseptic solutions.

Organization	Details
AORN, 2023	Skin antiseptics should be single-use containers (Page 651, 4.4.2) ¹
APIC	“APIC believes that preoperative skin preparations should be manufactured and marketed in single-use containers... we do not support use of multiple-use containers. Single-use containers should be designed in such a way that multiple-use is not possible.” ²
ORNAC, 2019	Single use containers of skin antiseptic shall be used for each patient. Single use containers decreases the risk of cross-contamination. (Page 2-41, 2.17.22) ³
FDA Drug Safety Communication, 2013	Recommends antiseptics for preoperative or pre-injection skin preparation be packaged in single-use containers. Applicators and any unused solution should be discarded after the single application. ⁴
INS, 2021	Use a single-use sterile applicator containing sterile solution, not a multiple use product (eg, bottle of antiseptic solution). (IV) Standard 33 ⁵

Solventum Surgical Skin Antisepsis Options

3M™ DuraPrep™ Surgical Solution

(Iodine Povacrylex [0.7% available iodine] and Isopropyl Alcohol, 74% w/w) Patient Preoperative Skin Preparation




3M™ SoluPrep™ S Sterile Antiseptic Solution

chlorhexidine gluconate (2% w/v) and isopropyl alcohol (70% v/v)
Patient Preoperative Skin Preparation



3M™ DuraPrep™ Surgical Solution Patient Preoperative Skin Preparation

- Dries to a water-insoluble film¹
- Resists removal by blood and irrigating solutions²
- Immobilizes bacteria remaining on the skin after prepping³
- Keep resident bacteria counts low for at least 48 hours* after blood and saline challenge⁴
- Applied in a single painted coat
- Provides statistically significant better drape adhesion for incise drapes than Chloraprep™ preoperative skin preparation⁵




*under simulated surgical conditions


1. Solventum data on file (LIMS 8198)
2. Solventum data on file(05-010565 and LIMS 8198)
3. Solventum data on file(05-004891)
4. Solventum data on file(05-010565)
5. Solventum data on file(05-010262)

Application Method

Beginning at the incision line and working your way out toward the periphery, paint a single uniform coat of solution on the skin.




One application method for dry and moist surgical sites




Note: Prior to the use of 3M™ DuraPrep™ surgical solution consult detailed directions contained in the Drug Facts.

Coverage Area



3M™ DuraPrep™ Surgical
Solution
26 mL
Covers 15" x 30" area
(450 sq inches)



BD ChloraPrep™ Patient
Preoperative Skin Preparation*
26 mL
Covers 13.2" x 13.2" area
(174 sq inches)

*BD Drug Facts accessed June 7, 2024
Actual percentage of torso coverage will vary by individual base on height and weight
BD and ChloraPrep are trademarks of Becton, Dickinson and Company

3M™ SoluPrep™ S Sterile Antiseptic Solution

chlorhexidine gluconate (2% w/v) and isopropyl alcohol (70% v/v)

Patient Preoperative Skin Preparation



Provides rapid, broad-spectrum antimicrobial activity* that persists for at least 96 hours**¹



FDA-approved sterile solution and sterile applicator^{2,3}



Visible on a variety of skin tones⁴



Contains an ingredient that can improve incise drape adherence^{†5}



Greater coverage area^{‡6}



Designed with workflow efficiency and convenience in mind⁷

1. Solventum Data on File. EM-05-013953, EM-05-014624.

2. Solventum Data on File. SV-EO-05-679401.

3. Solventum Data on File. SV-HEAT-05-764865.

4. Solventum Data on File. TECH-REPORT-05-943704, EM-05-014753, EM-05-014733.

5. Olson LK, Morse DJ, Paulson JE, Bernatchez SF. Evaluation of Incise Drape Lift Using 2% Chlorhexidine Gluconate/70% Isopropyl Alcohol Preoperative Skin Preparations in a Human Volunteer Knee Model. *Arthroplasty Today*. 2022 Mar 27;15:24-28. .

6. Solventum Data on File. EM-05-014815.

7. Solventum Data on File. DC-VER-05-726184, TECH-REPORT-05-742505, TECH-REPORT-05-757564.




*The clinical significance of in vitro data is unknown

**Study was performed on healthy volunteers under simulated surgical conditions

†Compared to skin prepped with Chloraprep® with 3M™ Ioban™ 2 Antimicrobial Incise Drapes in a simulated knee surgery model, 2019.

‡Compared to Chloraprep®, 26 mL and 10.5 mL, 2021.

Applicator Sizes and Coverage Area¹




Surgical Skin Preps and Incise Drape Adherence




OR



Incise Drapes





Surgical prep =
Reduce microbes

Surgical drapes =
Sterile field

Incise drape =
Sterile surface

3M™ Ioban™ 2 Antimicrobial Incise Drape



3M™ Ioban™ 2 Antimicrobial Incise Drape features



- Creates a sterile surface at the beginning of surgery




- Eliminates the need for towel clips and clamps to hold drapes in place



- Has 2% iodine that cannot be washed off it's in the adhesive



- Provides a barrier to liquids that can seep under the drape edges




- Film is breathable



- Helps prevent surgical skin preps from being washed away

Risk of contamination

There are always residual microbes that survive on the skin surface, in deeper skin layers, and in hair follicles



3M™ Ioban™ CHG Chlorhexidine Gluconate Incise Drape



Features & Benefits

CHG adhesive film^{*1}


Provides continuous, broad-spectrum antibacterial protection

Conformable and easy to use¹

Adheres to a variety of body contours and allows for limb manipulation.¹ It is easy to apply and remove

Trusted adherence

Provides adherence to skin comparable to Ioban 2 Antimicrobial Incise Drape¹



Breathable¹

Helps release moisture to encourage strong adherence

Withstands surgical challenges^{1**}

Such as irrigation fluid, limb manipulation, retraction, and in the presence of blood and fatty tissue. Protects against skin prep wash-off

Visibly transparent¹

Clear incise drape allows for easy viewing of the surgical site

1. Solventum Data on File, LAB-SUPPORT-05-867032
2. Solventum Data on File, SV-EO-05-867332, TECH-REPORT-05-889166

^{*}In vitro data, clinical significance is unknown.
^{**}In simulated surgery

Global guidelines: surgical incise drapes

ORNAC⁴

- “The use of adhesive surgical incise drapes without antimicrobial properties is not recommended to prevent SSIs”.
- “Iodophor impregnated adhesive drapes may be used in accordance with the manufacturer’s IFU to decrease SSIs, unless contraindicated by a patient’s allergy”.

KRINKO⁵

Evidence indicates antimicrobial-impregnated incise drapes result in reduction in bacterial colonization of the surgical site.

While bacterial colonization of the incision may predispose to subsequent SSI/PJIs, there is no literature to demonstrate that the use of incise drapes results in clinical differences in the rates of subsequent PJIs.

• APSIC⁶

“When using adhesive incise drapes, do not use non-iodophor-impregnated drapes for surgery as they may increase the risk of surgical site infection.”

“In orthopedic and cardiac surgical procedures where adhesive incise drapes are used, consider using an iodophor-impregnated incise drape, unless the patient has an iodine allergy or other contraindication.”

4. Operating Room Nurses Association of Canada. (2025, April). *Guidelines for perioperative Practice in Canada* (17th ed.). <http://www.ornac.ca/guidelines.phtml>

5. KRINKO Surgical Site Infection Prevention Guidelines, 2018.

6. APSIC Asia Pacific Society of Infection Control for The Prevention of Surgical Site Infections 2019.

Global guidelines: surgical incise drapes

NICE¹

“Do not use non-iodophor-impregnated incise drapes routinely for surgery as they may increase the risk of surgical site infection.”

“If an incise drape is required, use an iodophor-impregnated drape unless the patient has an iodine allergy.”

NHMRC²

“If an incise drape is required, use an iodophor-impregnated drape unless the patient has an iodine allergy.”

“Do not use non-iodophor-impregnated incise drapes routinely for surgery as they may increase the risk of surgical-site-infection.”

ACORN³

“Adhesive drapes with antimicrobial properties can be used in the critical aseptic field unless contraindicated (i.e. patient allergy).”

“These include but are not limited to iodophor-impregnated adhesive drapes.”

1. NICE National Institute for Health and Care Excellence . (2019, April 11). *Surgical site infections: prevention and treatment*. NICE National Institute for Health and Care Excellence . <https://www.nice.org.uk/guidance/ng125/chapter/Recommendations#intraoperative-phase>


2. NHMRC National Health and Medical Research Council guidelines for the prevention and control of infection in healthcare. NHMRC. 2019.

3. Australian College of Perioperative Nurses Ltd (ACORN) The New ACORN Standards. Volume-3-2023 Standards for safe and Quality Care in the Perioperative Environment (SSQCPE) for Organisations. Adelaide, South Australia: ACORN; 2023. Asepsis Standard, Critical Aseptic Field Maintenance.

Increasing the resistance of the patient

The SSI risk equation

**Risk
of
SSI** =
$$\frac{\text{Dose of bacteria} \times \text{Virulence of bacteria}}{\text{Resistance of the host (patient)}}$$



Control what you
can control



Mangram AJ, Horan TC, Pearson ML, Silver LC, and Jarvis WR. Guideline for prevention of surgical site infection. *Infect Control Hosp Epidemiol.* 1999;4:247-278. Retrieved from <https://www.ncbi.nlm.nih.gov/pubmed/24799638>.

The SSI risk equation

Risk of SSI

$$= \frac{\begin{array}{c} \boxed{5} \\ \text{Dose of bacteria} \end{array} \times \begin{array}{c} \boxed{5} \\ \text{Virulence of bacteria} \end{array}}{\begin{array}{c} \boxed{10} \\ \text{Resistance of the host (patient)} \end{array}} = 2.5$$

Control what you
can control



Mangram AJ, Horan TC, Pearson ML, Silver LC, and Jarvis WR. Guideline for prevention of surgical site infection. *Infect Control Hosp Epidemiol.* 1999;4:247-278. Retrieved from <https://www.ncbi.nlm.nih.gov/pubmed/24799638>.

Maintaining Normothermia

$$\text{Risk of SSI} = \frac{\text{Dose of bacteria} \times \text{Virulence of bacteria}}{\text{Resistance of the host (patient)}}$$


The diagram illustrates the risk of Surgical Site Infection (SSI) as a function of bacterial dose, virulence, and host resistance. The equation is presented as:

$$\text{Risk of SSI} = \frac{\text{Dose of bacteria} \times \text{Virulence of bacteria}}{\text{Resistance of the host (patient)}}$$

The numerator consists of two factors: "Dose of bacteria" (represented by an icon of a syringe and bacteria) and "Virulence of bacteria" (represented by an icon of various bacteria). The denominator is "Resistance of the host (patient)", represented by an icon of a patient in a hospital bed with a thermometer. This denominator is circled in red, indicating its critical role in maintaining normothermia and reducing SSI risk.

Guidelines:

Temperature management & Surgical site infection (SSI) prevention

World Health Organization¹ (WHO)

“Suggests the use of warming devices in the operating room for patient body warming with the purpose of reducing SSI.”

[Conditional recommendation, moderate quality of evidence]

Centers for Disease Control and prevention² (CDC)

“Maintain perioperative normothermia”

[Recommendation high to moderate-quality evidence]

Association of periOperative Registered Nurses³ (AORN)

“Implement preoperative warming (i.e., prewarming) as determined by the perioperative patient temperature management plan.” *[Recommendation]*

“Implement one or more active warming methods to prevent and treat inadvertent perioperative hypothermia as determined in the perioperative patient temperature management plan for all patients during all phases of perioperative care (ie, preoperative prewarming, intraoperative, postoperative).” *[Recommendation]*

“Implement intraoperative warming using methods defined in the perioperative temperature management plan.” *[Recommendation]*

Not a complete list of guidelines

1. Global guidelines for the prevention of surgical site infection, second edition. Geneva: World Health Organization; 2018. Licence: CC BY-NC-SA 3.0 IGO.
2. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention Guideline for the Prevention of Surgical Site Infection, 2017. *JAMA Surg.* 2017;152(8):784–791. doi:10.1001/jamasurg.2017.0904.
3. AORN. <https://aornguidelines.org/guidelines/content?sectionid=173731777&view=book>. Accessed April 25, 2025.

Guidelines:

Temperature management & Surgical site infection (SSI) prevention continued

American College of Surgeons¹

- “Studies have shown that intraoperative hypothermia is associated with increased risk of SSI, therefore, intraoperative maintenance is recommended.”
- “The use of preoperative warming before short clean cases has been shown to reduce SSI and is recommended.”
- “For longer cases, both preoperative warming and ongoing temperature monitoring and warming measures are recommended.”

1. Ban KA, Minei JP, Laronga C, Harbrecht BG, Jensen EH, Fry DE, Itani KMF, Dellinger EP, Ko CY, Duane TM. American College of Surgeons and Surgical Infection Society: Surgical Site Infection Guidelines, 2016 Update. J Am Coll Surg 2017;224:59-74.

Society for Healthcare Epidemiology of America/Infectious Diseases Society of America (SHEA/IDSA)²

- “Maintain normothermia (temperature of 35.5°C or more) during perioperative period.”
- “Randomized controlled trials have shown the benefits of both preoperative and intraoperative warming to reduce SSI rates and to reduce intraoperative blood loss”

[Recommendation quality of evidence grade I High]

2. Anderson DJ, Podgorny K, et al. Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. Infection Control and Hospital Epidemiology. 2014;35(6)

American Society of PeriAnesthesia Nurses (ASPAN)³

- “Identify patients at increased risk for hypothermia in the perioperative period.
[Recommendation]
- “Initiate preoperative/preprocedural nursing interventions to promote patient warming”
- “Patients with a temperature less than 36°C should have forced air warming initiated”
[Recommendation]
- “Refer to AORN’s Guideline for Prevention of Hypothermia for Intraoperative recommendations”

3. (2024). 2025-2026 Perianesthesia Nursing Standards, Practice Recommendations and Interpretive Statements (1st ed.). American Society of PeriAnesthesia Nurses.
<https://www.r2library.com/Resource/Title/9798895903001>

Not a complete list of guidelines

Copyright to all standard/guidelines owned by the respective association/society. No endorsement is implied or given.

Enhanced Recovery After Surgery: Guidelines

Association of periOperative Registered Nurses¹ (AORN)

- “Implement an ERAS program under the direction of the interdisciplinary team.” *[Recommendation]*
- **“Start patient warming in the preoperative period.”** *[Recommendation]*
- **“Implement measures to maintain intraoperative normothermia.”**
- *[Recommendation]*

1. *Implementation of ERAS (NEW)*. Aornguidelines.org. (n.d.).
<https://aornguidelines.org/guidelines/content?sectionid=246005266&view=book#246005266>

American Association of Nurse Anesthesiology² (AANA)

- **Preoperative**- “Maintain normothermia use warming gown”
- **Intraoperative** – “Maintain normothermia and use active warming to decrease oxygen demand from shivering post-operatively, as well as improve healing and decrease the risk of surgical site infection.”

2. *Enhanced Recovery After Surgery*. Aana.com. (n.d.).
<https://www.aana.com/practice/clinical-practice/clinical-practice-resources/enhanced-recovery-after-surgery/>

Not a complete list of guidelines

Copyright to all standard/guidelines owned by the respective association/society. No endorsement is implied or given.


Effects of anesthesia

Anesthesia drugs alter the hypothalamus' ability to maintain normothermia.

Anesthetic-induced vasodilation allows warmer blood from the core to flow freely and mix with the blood from the cooler periphery; cooler blood returns to the heart where it causes a drop in core temperature.^{1,2}

This is known as **heat redistribution**, commonly referred to as **RTD (redistribution temperature drop)**.

Anesthetized surgical patients cannot regulate temperature through behavioral changes and must rely on altered thermoregulation and clinical interventions to regulate temperature.^{1,2}



1. Sessler DI. Mild Perioperative Hypothermia. New Engl J Med. 1997;336(24):1730-1737.

2. Sessler DI. Temperature Monitoring. In: Miller RD, ed. Anesthesia. 3rd ed. New York: Churchill/Livingstone 1990

What is prewarming and why is it recommended?



The application of heat **prior to anesthesia** for the purpose of increasing total body temperature

Prewarming = “banking heat”

Why Is Prewarming Important?

It is difficult to **treat** RTD, but it can be **minimized** by prewarming.

Prewarming is the imposition of a temporary heat imbalance by active warming prior to anesthetic (general or neuraxial) induction to **increase the peripheral thermal compartment and mean body temperature without changing core body temperature**. Effective prewarming minimizes the decrease in core body temperature caused by the postinduction decrease of systemic vascular resistance (SVR) and the redistribution of body heat from the central to peripheral thermal compartment.¹ This is commonly referred to as Redistribution Temperature Drop (RTD).



1. Sessler DI, Schroeder M, Merrifield B, Matsukawa T, Cheng C. Optimal duration and temperature of prewarming. *Anesth.* 1995;82(3):674-681

3M™ Bair Hugger™ Warming Units

Model 775



- Hose-end temp sensing
- MERV 14 filter
- 2 fan speeds

Model 675



- Removable hose
- MERV 14 filter
- Filter change indicator
- Audible and visual over-temp and under-temp

Model 875




- Patient-adjustable controller
- Dust filter only
- Many mounting options

3M™ Bair Hugger™ Warming Gowns

3 Gown Styles – globally available (used with 875 warming unit)




Warming gowns debuted in 2003




OR gown ('01 Gown)

- Single comfort/ prewarming port (875)
- Single OR/clinical port (675/775) – Left=Large



Plus gown ('02 Gown)

- Single comfort/ prewarming port (875)
- Dual OR/clinical port (675/775) Left=Large



Flex gown ('03 Gown)

- Single comfort/ prewarming port (875)
- 2 Bair Hugger blankets inside (675/775)
 - Upper body insert with clinical port
 - Lower body insert with clinical port
- Left=Large

3M™ Bair Hugger™ Warming Gowns - Features



Hand inserts to warm hands



Easy side tie provides full coverage for modesty



Hook and Loop for easy opening at neck and shoulder



Secure hose attachment tabs on the hose end



Designed to be used through the entire perioperative journey



Patient-controlled temperature setting



Gown sides can be tucked during surgery to contain heat



Be sure to untie back ties before positioning patient on the OR table

3M™ Bair Hugger™ Full Access Underbody Blanket - 63500

- Placed and secured on OR table before the patient arrives in the room.
- Accommodates patients positioned supine, prone and lateral
- Can be used in other areas of the facility as well



Air is rerouted via quilted airflow channels when there is weight from the patient's body applied, to prevent warming on pressure points.



Always place a safety strap on surgical patients. Warming blankets are not designed to secure patients on OR table. Safety strap can be applied through the perforated portion on both sides of the blanket.



Full Access Underbody Blankets

- Pass-through slit design allows for flexible patient positioning and drawsheet use
- Fluid outlets minimize pooling of fluids on the surface of the blanket
- Soft, radiolucent materials accommodate imaging requirements
- Perforated area at head of blanket can be removed if desired (prone position)
- Adhesive strips and tuck flaps secure the blanket to the OR/procedure table
- Resealable hose ports at either end of the blanket provide options for hose placement
- One clear plastic head drape (included) helps retain warm air around the patient's head
- Outer edges of blanket draw up next to patient, surrounding them in warmth
- Micro-perforations at the foot end are closed to prevent direct heat application to the feet

3M™ Bair Hugger™ Spinal Underbody Blanket - 57501


Foot sling



Foot board




Pass the blanket to the underside of the patient support pads and back to the surface of the patient foot boards.



3M™ Bair Hugger™ Lithotomy Blanket - 58501

Features

- Underbody design provides full, unrestricted access to the patient and accommodate lateral, supine and lithotomy positions
- Fluid outlets minimize pooling of fluids on the surface of the blanket
- Pass-through slit design allows for flexible patient positioning and accommodates the use of a draw sheet
- Tape and tuck flaps secure the blanket to the table
- Generously-sized perineal cut-out at the base of the blanket provides clinician unobstructed access to patient
- Integrated tie strips can be loosely secured around the stirrups/leg supports
- The blanket draws up next to patients when inflated to surround them in warmth



Fluid drainage sites



Tuck flaps



Side perforations



Perineal perforated cut-out




Integrated tie strips



Always place a safety strap on surgical patients. Warming blankets are not designed to secure patients on OR table. Safety strap can be applied through the perforated portion on both sides of the blanket.

Temperature Modalities



AORN

- When possible, use the same modality throughout the surgical journey¹
- Zero-heat-flux can be used to measure core temperature¹

ASPAN

- Consider the ability to use the same modality in all phases of care²
- Zero-heat-flux (ZHF) thermometry provides an accurate reflection of the core temperature²

1. AORN. Recommended Practices for the Prevention of Unplanned Patient Hypothermia. AORN, Inc. 2023

2. ASPAN. Perianesthesia Nursing Standards, Practice Recommendations and Interpretative Statements. Practice Recommendation: Promotion of Normothermia in the Adult Patient. American Society of PeriAnesthesia Nurses. 2021-2022.

MHAUS

- Patients with general anesthesia >30 minutes should have core temperature monitored.¹

ACCM/IDSA

- Most accurate and reliable method should be chosen.²

1. MHAUS Recommendations. MHAUS Organization. <https://www.mhaus.org/healthcare-professionals/mhaus-recommendations>. Accessed June 5, 2023
2. O'Grady NP, Barie PS, Bartlett JG, Bleck T, Carroll K, Kalil AC, Linden P, Maki DG, Nierman D, Pasculle W, Masur H; American College of Critical Care Medicine; Infectious Diseases Society of America. Guidelines for evaluation of new fever in critically ill adult patients: 2008 update from the American College of Critical Care Medicine and the Infectious Diseases Society of America. Crit Care Med. 2008 Apr;36(4):1330-49.

NICE


- Use a modality that is a direct measurement or direct estimate of core temperature.¹
- Zero-heat-flux may be used to monitor core temperature.¹

ASA

- Continually monitor temperature when giving anesthetics.²

1. The Management of Inadvertent Perioperative Hypothermia in Adults. NICE. <https://www.nice.org.uk/guidance/cg65/chapter/Recommendations>. Accessed June 5, 2023
2. Standards for basic anesthetic monitoring. American Society of Anesthesiologists. <https://www.asahq.org/standards-and-guidelines/standards-for-basic-anesthetic-monitoring>. Amended 2015. Accessed June 5, 2023

Temperature Modality Accuracy



Wartzek T, Muhlsteff J, Imhoff M.
Temperature measurement. Biomedizinische
Technik/Biomedical Engineering.
2011;56(5):241-257


Overview



Non-invasive and accurate core temperature monitoring system that continuously measures temperature with a single-use sensor.




How It Works



This image shows core temperature rising to the surface through an isothermal pathway.

The Sensor Unveiled



Above image shows the inverse direction of how the product appears when in-use.



Control Unit



Current temperature is displayed in large numbers

Previous temperature data is displayed as a trend graph

Temperature trend graph displays previous two hours of data

Temperature can be displayed in Celsius or Fahrenheit

Data Transfer

Temperature data may be continuously transmitted to a compatible patient vital signs monitor and electronic medical record.



3M™ Bair Hugger™
Temperature Monitoring
System



Patient Vital Signs Monitor



EMR

3M™ Bair Hugger™ Temperature Monitoring System

Non-invasive and accurate core temperature monitoring system that continuously measures temperature with a single-use sensor.



Rx Only. This product is intended to be used by trained medical professionals in a clinical/surgical setting.



Thank you!

 solventum